Cleaning and Care Instructions	Instrucciones de cuidado y limpieza
Anleitung zur Reinigung und Pflege	Instructions d'entretien et de nettoyage
Istruzioni per la cura e pulizia	Instruções de limpeza e cuidado
Návod k čištění a údržbě	Οδηγίες καθαρισμού και φροντίδας
Reinigings- en onderhoudsinstructies	Instrukcje czyszczenia i pielęgnacji
Temizleme ve Bakım Talimatları	Instruktioner för rengöring och skötsel

Montgomery® Safe-T-Tube™

REF: _____

LOT: _____

IMPORTANT:

Please deliver these important instructions to the Safe-T-Tube patient / caregiver.

English	3
Español	12
Deutsch	22
Français	32
Italiano.....	42
Português	52
Český	62
ελληνικά.....	71
Nederlands	82
Polski.....	92
Türkçe	102
Svenska.....	111

1 Introduction

1.1 Dear Patient,

The tracheostomy tube you have in your throat is a *Montgomery Safe-T-Tube* (also: Tracheal T-Tube or short: T-Tube). Unlike standard tracheostomy tubes, the Montgomery Safe-T-Tube does not have an inner cannula to be removed for cleaning.

T-Tubes have been used successfully to treat patients like you for many years.

However, care must be taken to keep the T-Tube clean and clear of mucus / secretions which, if left unattended, could eventually result in partial or complete blockage of the T-Tube.

Therefore, please read and follow these Cleaning and Care Instructions and save these instructions. If after reading these instructions you have any questions or concerns, please contact your physician's office, Boston Medical Products or your BMPI Distributor.

1.2 Symbols Glossary

Symbol	Description
	Caution: Consult Instructions for Use
	Do not use if package is damaged
	Keep away from direct sunlight
	Keep dry
	Use-by date
	Catalog number
	Batch code
	Quantity per packaging unit

Symbol	Description
	Manufacturer
	(EU) Authorized representative in the European Community
	(USA) Caution: Federal Law restricts this device to sale by or on the order of a physician.

Table 1: Symbols Glossary

1.3 Safety Information Marking

WARNING

Non-compliance may result in serious injuries, serious deterioration of your general condition or your death.

CAUTION

Non-compliance may result in light or moderate injuries or a light or moderate deterioration of your general condition.

2 Important Safety Information

WARNING

- Read these instructions carefully, adhere to them and save them. Otherwise there are risks to your health.
- Always keep the T-tube clean and free from obstructions. Otherwise there is a risk for you to suffocate.
- To be prepared for an emergency removal of your T-Tube: Always keep a standard tracheostomy tube in a size smaller than that of your T-Tube at hand. Otherwise there is a risk for you to suffocate.

3 Your T-Tube

3.1 Components

These items are components of your T-Tube and need to be handed to you from your physician:

- 1 x Patient Pak [™]
 - 1 x Cleaning and Care Instructions (N-STCC)

- 1 x Plug/Ring Set storage container
- 2 x Plug/Ring-Set, of which 1 mounted on the T-Tube
- 1 x implant card

3.2 Replacement Parts

Available replacement parts:

- Plug/Ring-Set

Please contact Boston Medical Products or your BMPI Distributor.

4 Prior to Leaving the Clinic

The production of mucus and secretions varies between the patients. Therefore the required cleaning techniques and their frequency differ individually. Also, the condition of the skin around the tracheostoma differs from patient to patient.

Therefore please check with your physician prior to leaving the clinic:

- How often do you need to apply these cleaning techniques?
- Where do you obtain the material required for the cleaning techniques?
- How and how often do you have to take care for the skin around your tracheostoma?

5 Configuration of the T-Tube

WARNING

- If the upper limb of your T-tube is blocked (e.g., by a foreign object or by concomitant use of a closed laryngeal stent): Ensure the external limb is free of obstruction and that you can safely breathe through this airway. Otherwise there is a risk for you to suffocate.

Provided the upper limb is free of obstruction and you can safely breathe through this airway: Occlude the external limb with the plug. If the external limb is not closed for a longer period of time: Supply humidified air.


Illustration 1: Left: Upper limb is free, Right: Upper limb is blocked

ATTENTION: As a general rule, the T-tube must be secured with the safety ring. If the T-tube cannot be secured with the safety ring for a short time (e.g. during tracheostoma care): Hold the external limb firmly so that the T-tube does not slip into the trachea.

6 Cleaning And Care of Your T-Tube

6.1 Required Equipment and Materials

- Standard tracheostomy tube with obturator, smaller than the T-Tube used
- Isotonic saline solution
- Pipette
- Flexible tracheostomy tube brush with a bristled, nonabrasive distal tip
- Suction unit
- Suction catheter, size: 8 - 12 Ch
- Humidification equipment (if recommended by the physician)
- Hemostat

6.2 Lavage

⚠ CAUTION

- Do not inhale while the saline solution is inside your T-Tube. Otherwise there is the risk that you inhale the saline solution and suffer from shortness of breath.

During lavage you squirt saline solution into your T-Tube and then cough to clean the inside of your T-Tube.

Sometimes suctioning is also required to remove excess mucus. [▶
Suctioning, page 7]

1. Take 5 - 10 ml of saline solution with the pipette.
2. Squirt the saline solution into your T-Tube.
CAUTION: Do not inhale while the saline solution is inside your T-Tube.
3. Close your T-Tube with a finger and cough.
4. Repeat if necessary.

TIP: If your pipette does not have a scale, you can use a teaspoon to determine the amount of saline solution you need: Fill your pipette with the amount of saline solution equivalent to 1 - 2 teaspoons and note the fill level on your pipette.


6.3 Suctioning

CAUTION


- Do not insert the suction catheter beyond the ends of your T-Tube.

Otherwise you might suffer from irritations / coughing.

1. Insert the suction catheter gently into the outer limb of your T-Tube.

2. Gently bend the external limb of your T-Tube up to direct the suction catheter down.
3. Suction the lower limb of your T-Tube.
4. Remove the suction catheter slightly.

5. Gently bend the external limb of your T-Tube down to direct the suction catheter up.
6. Suction the upper limb of your T-Tube.
7. Remove the suction catheter.

TIP: On the tube of the suction catheter, mark how far you have to insert the suction catheter into the upper / lower limb of your T-Tube.

7 Emergency Measures

7.1 Signs for a Blockage of the T-Tube / for Troubled Breathing

- Noise from the T-Tube
- Difficulty of speaking
- Anxiety
- Crowing sounds from the throat
- Change in skin color to pale or blue
- Skin feels cold, but damp to touch

If the patient shows any signs of troubled breathing, act immediately:

7.2 Assessing the Patient's Situation

1. Place the palm of your hand over the T-Tube to assess air movement.
2. Listen for equal breath sounds.
3. If troubled breathing is confirmed or if no air movement can be confirmed: Make an emergency call.
4. After that act depending on the situation:
 - [▶ If Troubled Breathing is Confirmed, page 8]
 - [▶ If No Breathing Can Be Confirmed, page 9]

7.3 If Troubled Breathing is Confirmed

1. Suction the T-Tube. [▶ Suctioning, page 7]

2. Re-assess the patient's breathing.
3. Call the physician to alert him/her of the situation.
4. Upon arrival of the emergency personnel: Inform the emergency personnel that there is a T-Tube in place and show them these instructions.

ATTENTION: In the waiting time, stay with the patient and continue to perform suctioning.

7.4 If No Breathing Can Be Confirmed

WARNING

- Do not grasp the T-Tube with pointed or sharp instruments to remove it.
Otherwise, the T-tube may be damaged and parts of the T-tube may enter the lower airways. There is a risk of suffocation.

If no breathing can be confirmed the T-Tube might be completely blocked. In this case the T-Tube must be removed immediately and be replaced with a a smaller size standard tracheostomy tube:

1. Prepare the standard tracheostomy tube and a hemostat.
2. With one hand, grasp the external limb of the T-Tube. At the same time hold the skin above and below the T-Tube with two fingers of the other hand.
3. Remove the T-Tube with a steady, slow pull. If the external limb of the T-Tube tears off: Grasp the remaining portion of the T-Tube with the hemostat and remove it.
4. Insert the standard tracheostomy tube and remove the obturator.
5. If necessary: Suction the airway and ventilate the patient through the standard tracheostomy tube.
6. Upon arrival of the emergency personnel: Inform the emergency personnel that there was a T-Tube in place and show them these instructions.

8 Warranty

The reliability of the product's material and design at the time of shipment is guaranteed. The manufacturer does not know either the diagnosis of the patient or the nature of the application and has no influence on the conditions under which the product is used. The storage conditions after delivery of the product are also beyond the manufacturer's area of responsibility.

Due to biological and individual differences, no product is 100% effective under all circumstances.

Therefore, the manufacturer cannot guarantee a positive effect or the absence of negative effects for product application. The medical staff must use the product on the basis of their medical training and experience, and they are responsible for correct application.

The warranty (repair or replacement) applies only if the product is used in accordance with these Instructions for Use (for instruments, particularly with regard to handling, cleaning, sterilization and maintenance); the warranty period starts on the delivery date.

If you have reason to believe that a new product is faulty, please contact the Customer Service in writing immediately and provide as detailed a description as possible of the fault, the REF (product code), and the LOT (batch code) and/or series number. All allegedly defective products must be returned to us for inspection. Instruments have to be completely cleaned and sterilized, appropriate documentation must be enclosed with the return.

If the manufacturer finds that despite all due care the product was defective at the time of delivery, he will repair the product or replace it promptly. If repair or replacement of the product is not possible, the buyer has the right to cancel the purchase or to reduce the payment, but by a maximum of the purchase price amount.

Additional claims or those not mentioned here due to defect, and other claims regardless of the legal reason, including those based on illegal acts and for compensation of immaterial damages against the manufacturer, his agents, dealers and suppliers, are excluded unless existing law is contrary to the liability exclusion, e.g. in cases of intent or gross negligence or in the event of physical injury.

All claims based on the consequences of non-compliance with the Instructions for Use, including specified indications, contraindications, warnings, instructions, application, storage and off-label use, as well as the consequences of a combination with third-party products are excluded.

Furthermore, all claims that result from the use of products that have expired, or were used despite the obvious damage to the packaging, or resterilized and/or recycled contrary to the Instructions for Use, are excluded.

No one is allowed to change the above conditions, make further warranty or liability declarations, or guarantee any properties that surpass those specified in the Instructions.

The General Terms and Conditions of the manufacturer, which can be accessed at <http://www.bosmed.com> apply in all remaining instances.

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Boston Medical Products, Inc.