

Diameters and lengths - straight stents

DUMON™ BD / DUMON™ TD / DUMON™ TF

Determine the desired stent length (stenosis length + 10 mm) and diameter and select from the available stents.



Custom-made products are possible.
Contact us!

Length (mm) ►	20	30	40	50	60	70	80	90	100	110
Outer Diameter (mm) ▼										
9	02841S(BD)	04341S(BD)	04371S(BD)	04391S(BD)	04411S(BD)	04421S(BD)				
10	03721S(BD)	03741S(BD)	03761S(BD)	03781S(BD)	03801S(BD)	03821S(BD)				
11	03001S(TD) 03841S(BD)	03021S(TD) 03861S(BD)	03041S(TD) 03881S(BD)	03061S(TD) 03901S(BD)	03081S(TD) 03921S(BD)	03101S(TD)				
12	03121S(TD) 03961S(BD)	03141S(TD) 03141FS(TF)* 03981S(BD)	03161S(TD) 03161FS(TF)* 04001S(BD)	03181S(TD) 04021S(BD)	03201S(TD) 04041S(BD)	03221S(TD) 04051S(BD)	03231S(TD)			
13		03261S(TD) 03261FS(TF)	03281S(TD) 03281FS(TF)	03301S(TD) 03301FS(TF)	03321S(TD) 03321FS(TF)	03341S(TD) 03341FS(TF)				
14		03381S(TD) 03381FS(TF)	03401S(TD) 03401FS(TF)	03421S(TD) 03421FS(TF)	03441S(TD) 03441FS(TF)	03461S(TD) 03461FS(TF)				
15		03531S(TD)	03541S(TD)	03551S(TD)	03561S(TD)	03571S(TD)	03581S(TD)	03591S(TD)	03491S(TD)	03501S(TD)
16		03621S(TD)	03641S(TD) 03641FS(TF)	03661S(TD) 03661FS(TF)	03681S(TD) 03681FS(TF)	03701S(TD) 03701FS(TF)	04281S(TD) 04281FS(TF)	04291S(TD) 04291FS(TF)	04301S(TD) 04301FS(TF)	04311S(TD)
18			04141S(TD) 04141FS(TF)*	04151S(TD) 04151FS(TF)*	04161FS(TF)*	04171FS(TF)*	04181FS(TF)*	04191FS(TF)*	04201FS(TF)*	
20			04231FS(TF)*	04241FS(TF)*	04251FS(TF)*	04261FS(TF)*	04271FS(TF)*	04431FS(TF)*	04441FS(TF)*	04451FS(TF)*



DUMON™ BD



DUMON™ TD



DUMON™ TF

REF + R = Radiopaque *transparent only

DUMON™ stents are supplied sterile in a double sterile pouch, one stent per carton.

