


## Diameters and lengths - straight stents

### DUMON™ BD / DUMON™ TD / DUMON™ TF

Determine the desired stent length (stenosis length + 10 mm) and diameter and select from the available stents.

 Custom-made products are possible. Contact us!

| Length (mm) ▶         | 20                     | 30                                    | 40                                    | 50                       | 60                      | 70                      | 80                      | 90                      | 100                     | 110         |
|-----------------------|------------------------|---------------------------------------|---------------------------------------|--------------------------|-------------------------|-------------------------|-------------------------|-------------------------|-------------------------|-------------|
| Outer Diameter (mm) ▼ |                        |                                       |                                       |                          |                         |                         |                         |                         |                         |             |
| 9                     | 02841(BD)              | 04341(BD)                             | 04371(BD)                             | 04391(BD)                | 04411(BD)               | 04421(BD)               |                         |                         |                         |             |
| 10                    | 03721(BD)              | 03741(BD)                             | 03761(BD)                             | 03781(BD)                | 03801(BD)               | 03821(BD)               |                         |                         |                         |             |
| 11                    | 03001(TD)<br>03841(BD) | 03021(TD)<br>03861(BD)                | 03041(TD)<br>03881(BD)                | 03061(TD)<br>03901(BD)   | 03081(TD)<br>03921(BD)  | 03101(TD)               |                         |                         |                         |             |
| 12                    | 03121(TD)<br>03961(BD) | 03141(TD)<br>03141F(TF)*<br>03981(BD) | 03161(TD)<br>03161F(TF)*<br>04001(BD) | 03181(TD)<br>04021(BD)   | 03201(TD)<br>04041(BD)  | 03221(TD)<br>04051(BD)  | 03231(TD)               |                         |                         |             |
| 13                    |                        | 03261(TD)<br>03261F(TF)               | 03281(TD)<br>03281F(TF)               | 03301(TD)<br>03301F(TF)  | 03321(TD)<br>03321F(TF) | 03341(TD)<br>03341F(TF) |                         |                         |                         |             |
| 14                    |                        | 03381(TD)<br>03381F(TF)               | 03401(TD)<br>03401F(TF)               | 03421(TD)<br>03421F(TF)  | 03441(TD)<br>03441F(TF) | 03461(TD)<br>03461F(TF) |                         |                         |                         |             |
| 15                    |                        | 03531(TD)                             | 03541(TD)                             | 03551(TD)                | 03561(TD)               | 03571(TD)               | 03581(TD)               | 03591(TD)               | 03491(TD)               | 03501(TD)   |
| 16                    |                        | 03621(TD)                             | 03641(TD)<br>03641F(TF)               | 03661(TD)<br>03661F(TF)  | 03681(TD)<br>03681F(TF) | 03701(TD)<br>03701F(TF) | 04281(TD)<br>04281F(TF) | 04291(TD)<br>04291F(TF) | 04301(TD)<br>04301F(TF) | 04311(TD)   |
| 18                    |                        |                                       | 04141(TD)<br>04141F(TF)*              | 04151(TD)<br>04151F(TF)* | 04161F(TF)*             | 04171F(TF)*             | 04181F(TF)*             | 04191F(TF)*             | 04201F(TF)*             |             |
| 20                    |                        |                                       | 04231F(TF)*                           | 04241F(TF)*              | 04251F(TF)*             | 04261F(TF)*             | 04271F(TF)*             | 04431F(TF)*             | 04441F(TF)*             | 04451F(TF)* |



DUMON™ BD



DUMON™ TD



DUMON™ TF

REF + R = Radiopaque \*transparent only  
DUMON™ stents are supplied non-sterile, one stent per carton.

